


Accessory kit: Bottle OMNIS IC, PP, 2 L

6.05000.010


OMNIS IC accessory set to install the reagent bottle onto the Reagent Organizer / Reagent Organizer Liquid Control.

Phạm vi giao hàng 6.05000.010


Order Number	Description	Quantity
6.01602.100	Flaschenaufsatz GL 45, 2 x M6, UNF 10/32	1 PCS



6.01608.000	Bottle, plastic (PP), 2 L, GL 45	1 PCS
-------------	----------------------------------	-------




6.01834.010	FEP-Schlauch M6 mit Filter, 0.27 m	1 PCS
-------------	------------------------------------	-------



6.02050.000

Flaschenhalter, RO


1 PCS



6.02050.010

Flaschenbefestigung RO


1 PCS



Phụ kiện tùy chọn

6.05801.040

**Volumensensor, ROLC-
Kit**



6.02061.000

Reagenzorganisierer

