

Zubehör-Kit: OMNIS IC AS COOL

6.05000.060

Phạm vi giao hàng 6.05000.060


Order Number	Description	Quantity
6.01816.010	Silicone tubing, 8/12 mm, 2 m Drainage tubing	2 PCS



6.02041.510

OMNIS Probenrack 40 x 5 mL


2 PCS



6.02061.120

Tropfensammler, Extension


1 PCS



6.02743.027

Vial PS, 5 mL, 200 Stk.

1 PCS



6.02743.030

Bottle, PP, 50 mL, 12 pieces

1 PCS




The 50 mL bottles are made of polypropylene (PP) plastic and can be used for samples, standards, check standards and rinsing solutions. Production and packaging both take place under sterile and clean room conditions. This guarantees minimum contamination. The sample vessels are subjected to random checks for leachable anions and cations.

6.02743.087

Kappe PP, AD 11.5 mm, 200 Stk.

1 PCS



6.02040.010

OMNIS IC AS COOL Tray

1 PCS

