



NH₃-selective gas membrane electrode (low conc.)

6.0506.100

Combined ammonia-selective gas membrane electrode for low concentrations 10^{-2} mol/L, e.g.:

- in drinking water
- in boiler feed water
- in mineral water

If necessary, the complete membrane module can be easily exchanged. They are individually tested and certified.

Phạm vi giao hàng 6.0506.100

Qt.	Order no.	Description
-----	-----------	-------------

100 Stk

6.1236.050

Sleeve with SGJ 14/12 mm

Sleeve with SGJ 14/12 mm, polyethylene.



100 Stk

6.2008.040

Storage vessel

Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.



Phụ kiện tùy chọn

Order no.	Description	
6.1255.000	Replacement Kit for NH ₃ -selective Electrode 6.0506.100 Consisting of: replacement module (3 pieces), measuring electrolyte (50 mL)	
6.2104.020	Electrode cable / 1 m / F For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F).	
6.2316.030	Inner filling solution to 6.0506.100 electrode Measuring electrolyte for NH ₃ -selective electrode (50 mL)	