



Combined polymer membrane electrode, K

6.0510.110

Combined potassium-selective electrode with polymer membrane.
This ISE is suitable for:

- ion measurements of K^+ (10^{-7} to 1 mol/L) in aqueous solutions
- titrations in aqueous solutions

Thanks to a robust/break-proof plastic shaft made of propylene and an impact protection for the polymer membrane, this sensor is mechanically very resistant.

The reference electrolyte used is $c(CH_3COOLi) = 1$ mol/L.

Phạm vi giao hàng 6.0510.110

Qt.	Order no.	Description
-----	-----------	-------------

10 Stk

6.1236.050

Sleeve with SGJ 14/12 mm

Sleeve with SGJ 14/12 mm, polyethylene.



10 Stk

6.2008.040

Storage vessel

Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.



Phụ kiện tùy chọn

Order no.	Description	
6.2328.000	Lithium acetate electrolyte, 1 mol/L (50 mL) Lithium acetate electrolyte, 1 mol/L (reference electrolyte for combined K ⁺ -ISE 6.0510.110)	
6.2104.020	Electrode cable / 1 m / F For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F).	