



## OMNIS IC-Ausr.: PP Reg.


## Phạm vi giao hàng

| Order Number                                                                        | Description                                                                                                                                                               | Quantity |
|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|
| 6.02837.040                                                                         | Filterrohr                                                                                                                                                                | 1 PCS    |
|    |                                                                                                                                                                           |          |
| 6.05000.010                                                                         | <b>Accessory kit: Bottle OMNIS IC, PP, 2 L</b><br><br>OMNIS IC accessory set to install the reagent bottle onto the Reagent Organizer / Reagent Organizer Liquid Control. | 1 PCS    |
|  |                                                                                                                                                                           |          |
| 6.2744.034                                                                          | <b>Coupling olive/UNF 10/32, 2x</b><br><br>Connection between pressure screw and pump tubing. 2 pieces. For IC instruments with peristaltic pump.                         | 1 PCS    |
|  |                                                                                                                                                                           |          |

6.01803.110

**PTFE-Kapillare RO - PP, 900/0.5 mm**

1 PCS



---

6.01826.010

**Pumpschlauch, ID 0.5 mm, orange - gelb**

1 PCS

