




OMNIS IC-Ausr.: PCR, DD

Phạm vi giao hàng

Order Number	Description	Quantity
6.02744.010	flanschlose Verschraubung M6, 1/8", 2 Stk.	2 PCS
		
6.02837.040	Filterrohr	2 PCS
		
6.05000.010	Accessory kit: Bottle OMNIS IC, PP, 2 L OMNIS IC accessory set to install the reagent bottle onto the Reagent Organizer / Reagent Organizer Liquid Control.	2 PCS
		

6.2744.070

Pressure screw short

1 PCS


Short version. With UNF 10/32 connection. 5 pieces. For the connection of PEEK capillaries



6.02027.030

Clip Befestigungsschiene, 2 x AD 1/16

2 PCS



6.2744.150

T connector UNF 10/32

1 PCS

T connector for 3 capillaries.



6.1831.030

PEEK capillary 0.75 mm i.d. / 3 m

1 PCS


For all IC components.



6.01805.400

FEP-Schlauch M6, 1/8" x 1.59 mm, 3 m


1 PCS



6.02006.120

Halter, Befestigungsschiene


1 PCS



6.02061.130

Tropfensammler, Unit

1 PCS



6.02845.100

PCR-Kartusche, PTFE, ID 0.5 mm

1 PCS

