



OMNIS IC Ausr.: EPOD

Phạm vi giao hàng

Order Number	Description	Quantity
6.1621.000	Container 10 L As rinsing or waste container in automated systems.	1 PCS
		
6.02744.010	flanschlose Verschraubung M6, 1/8", 2 Stk.	3 PCS
		
6.05000.010	Accessory kit: Bottle OMNIS IC, PP, 2 L OMNIS IC accessory set to install the reagent bottle onto the Reagent Organizer / Reagent Organizer Liquid Control.	1 PCS
		

6.1805.530

FEP tubing / M6 / 2 m

1 PCS

With light and kink protection



6.2744.200

Adapter UNF 10/32 outer / M6 inner

1 PCS


Connecting 1/16 in. PEEK capillaries and transfer tubing (M6)



6.01602.110

Flaschenaufsatz S 40, 2x M6, UNF 10/32

1 PCS



6.01503.250

OMNIS 50 mL cylinder unit without accessories

1 PCS

Intelligent 50 mL cylinder unit for an OMNIS Titrator, Titration Module or Dosing Module. Without dosing tubing and antidiffusion tip.



6.01805.400

FEP-Schlauch M6, 1/8" x 1.59 mm, 3 m

2 PCS

