



Ag/AgCl reference system

6.0724.140

Silver / silver chloride reference system with SGJ 14/15 for the modular reference system in combination with the reference electrolyte vessel made of PTFE (6.1240.000) or glass (6.1240.020).

Phụ kiện tùy chọn

Order no.	Description
6.1240.000	Electrolyte vessel / diameter 3 mm With PTFE capillary and ceramic diaphragm (3 mm)



6.1240.020 Electrolyte vessel, diameter 1 mm
With ceramic diaphragm (D type).



6.2308.020 Electrolyte 3 mol/L KCl (250 mL)
Electrolyte solution $c(\text{KCl}) = 3 \text{ mol/L}$, (for Ag/AgCl reference systems)

