



Complete LL-Ag/AgCl(gel) reference electrode

6.0730.100

Maintenance-free Double-Junction Ag/AgCl reference electrode for voltammetry. Made of glass with electrolyte $c(\text{KCl}) = 3 \text{ mol/L}$ (gel), for pH 0 - 8. Suitable for use in continuous operation in automated systems and online instruments.

Phạm vi giao hàng 6.0730.100

Qt.	Order no.	Description
1 PCS	6.1236.050	Sleeve with SGJ 14/12 mm
		Sleeve with SGJ 14/12 mm, polyethylene.



Together with 6.2008.050 storage vessel holder, it provides a support for the electrode on 807 Dosing Units.



Phụ kiện tùy chọn

Order no.	Description
6.2112.130	Adapter plug contact 5.5 mm / screw contact, gold-plated For VA electrodes

