



Ag/AgCl reference electrode (length 18.5 cm)

6.0736.110

Silver / silver chloride reference electrode, installation length 16.8 cm. This reference electrode is well suited for automated applications. The flexible ground-joint diaphragm, which is insensitive to contamination, can be replaced at any time. The recommended reference electrolyte is $c(\text{KCl}) = 3 \text{ mol/L}$. This sensor is supplied without electrolyte filling.

Phụ kiện tùy chọn

Order no.	Description
6.2106.020	Strand / 1 m / 2 x B Electrode cable with 2 x plug B (4 mm) for connecting electrodes with plug-in head B



6.2106.060

Strand / 2 m / 2 x B

Electrode cable with 2 x plug B (4 mm) for connecting electrodes with plug-in head B



6.2308.020

Electrolyte 3 mol/L KCl (250 mL)

Electrolyte solution $c(\text{KCl}) = 3 \text{ mol/L}$, (for Ag/AgCl reference systems)



6.2310.010

Electrolyte KNO_3 1 mol/L 250 mL

Electrolyte solution KNO_3 1 M (reference electrolyte for combined silver electrodes)



6.2320.000

TEABr 0.4 mol/L in ethylene glycol (250 mL)

Electrolyte solution TEABr (tetraethylammonium bromide in ethylene glycol), $c(\text{TEABr}) = 0.4 \text{ mol/L}$

