



High Throughput Lab-grade Raman Probe (532 nm), 5m


607510090

High-throughput lab-grade fiber optic Raman probe assembly with trigger button. 5 meter fiber length.

Compatible with i-Raman NxG 532.

Phạm vi giao hàng 607510090

Qt.	Order no.	Description
1 PCS	607510430	Distance Regulator Cap
Distance regulator cap for standard lab-grade Raman probe shaft.		



Replacement shaft for 532 nm excitation lab-grade Raman probe assembly. Working distance 5.4 mm.

