



High Throughput Lab-grade Raman Probe (785 nm)

607510140


High-throughput lab-grade fiber optic Raman probe assembly with trigger button. 1.5 meter fiber length.

Compatible with i-Raman NxG 785S and 785H models.

Phạm vi giao hàng 607510140

Qt.	Order no.	Description
1 PCS	607510430	Distance Regulator Cap

Distance regulator cap for standard lab-grade Raman probe shaft.



Replacement shaft for 785 nm excitation lab-grade Raman probe assembly. Working distance 5.5 mm.

