



High Throughput Lab-grade Raman Probe (1064 nm), 5 m


607510270

High-throughput lab-grade fiber optic Raman probe assembly with trigger button. 5 meter fiber length.

Compatible with i-Raman NxG 1064.

Phạm vi giao hàng 607510270

Qt.	Order no.	Description
1 PCS	607510430	Distance Regulator Cap
Distance regulator cap for standard lab-grade Raman probe shaft.		



Replacement shaft for 1064 nm excitation lab-grade Raman probe assembly. Working distance 5.9 mm.

