



# High Throughput Lab-grade Raman Probe for i-Raman Duo


607510360

High-throughput lab-grade fiber optic Raman probe assembly with trigger button. 1.5 meter fiber length. Compatible with i-Raman Duo.

Phạm vi giao hàng 607510360

| Qt.   | Order no. | Description            |
|-------|-----------|------------------------|
| 1 PCS | 607510430 | Distance Regulator Cap |

Distance regulator cap for standard lab-grade Raman probe shaft.



Replacement shaft for 1064 nm excitation lab-grade Raman probe assembly. Working distance 5.9 mm.

