



High Throughput Lab-grade Raman Probe for i-Raman Duo, 1.5 m


6.07510.370

High-throughput lab-grade fiber optic Raman probe assembly with trigger button. 1.5 meter fiber length. Compatible with i-Raman Duo.

Phạm vi giao hàng 6.07510.370

Qt.	Order no.	Description
1 PCS	607510430	Distance Regulator Cap

Distance regulator cap for standard lab-grade Raman probe shaft.



Replacement shaft for 1064 nm excitation lab-grade Raman probe assembly. Working distance 5.9 mm.

