



Lab-Grade Raman Probe (785 nm), 5 m


607510460

Lab-grade fiber optic Raman probe assembly with trigger button. 5 meter fiber length. Compatible with i-Raman Plus 785S and 785HR models.

Phạm vi giao hàng 607510460

Qt.	Order no.	Description
1 PCS	607510430	Distance Regulator Cap

Distance regulator cap for standard lab-grade Raman probe shaft.



Replacement shaft for 785 nm excitation lab-grade Raman probe assembly. Working distance 5.5 mm.

