



Conductivity measuring cell $c = 0.9 \text{ cm}^{-1}$

6.0910.120

Conductivity measuring cell with cell constant $c = 0.9 \text{ cm}^{-1}$, without integrated temperature sensor. This sensor is suitable for measurements of medium conductivities (1 $\mu\text{S}/\text{cm}$ to 100 mS/cm), e.g., in:

- drinking water
- surface water
- wastewater

This sensor is suitable for automated applications (includes an SGJ 14/15).

Phụ kiện tùy chọn

Order no.	Description
-----------	-------------

6.2301.060 KCl ion standard 250 mL

KCl ion standard ($c(\text{ion}) = 0.1000 \pm 0.0005 \text{ mol/L}$), can also be used as standard solution for conductometry 12.87 mS/cm (25°C).



6.2324.010 Conductivity standard $100 \mu\text{S/cm}$ 250 mL

Conductivity standard for calibration of conductivity measuring cells with cell constant = $0.1/\text{cm}$.



6.2324.110 Conductivity standard $100 \mu\text{S/cm}$, 5 x 30 mL

Conductivity standard for calibration of conductivity measuring cells with cell constant = $0.1/\text{cm}$.



6.2104.020 Electrode cable / 1 m / F

For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F).



6.2104.030

Electrode Cable 2 m / F

For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F).

