




Conductivity measuring cell $c = 0.8 \text{ cm}^{-1}$ with Pt1000 (fixed cable)

6.0912.110

Conductivity measuring cell with cell constant $c = 0.8 \text{ cm}^{-1}$, with integrated Pt1000 temperature sensor and with fixed cable (4 x B (4 mm)). Thanks to the robust/break-proof plastic shaft made of polypropylene and an impact protection, this sensor is mechanically very resistant. This sensor is suitable for measurements of medium conductivities (1 $\mu\text{S}/\text{cm}$ to 100 mS/cm), e.g., in:

- drinking water
- surface water
- wastewater

The sensor can be connected to an 856 Conductivity Module with the adapter box (6.2103.160).

Phụ kiện tùy chọn

Order no.	Description
-----------	-------------

6.2301.060

KCl ion standard 250 mL

KCl ion standard ($c(\text{ion}) = 0.1000 \pm 0.0005 \text{ mol/L}$), can also be used as standard solution for conductometry 12.87 mS/cm (25°C).



6.2324.110

Conductivity standard $100 \mu\text{S/cm}$, 5 x 30 mL

Conductivity standard for calibration of conductivity measuring cells with cell constant = $0.1/\text{cm}$.



6.2324.010

Conductivity standard $100 \mu\text{S/cm}$ 250 mL

Conductivity standard for calibration of conductivity measuring cells with cell constant = $0.1/\text{cm}$.



6.2103.160

Adapter 4 x socket B – plug N

Adapter box for the connection of classical Metrohm Conductivity measuring cells with 4 banana plugs to the 856 Conductivity Module.

