




Conductivity measuring cell $c = 0.1 \text{ cm}^{-1}$ with Pt1000 (fixed cable)

6.0914.040

Conductivity measuring cell made of stainless steel with cell constant $c = 0.1 \text{ cm}^{-1}$, with integrated Pt1000 temperature sensor and with fixed cable (4 x B (4 mm)).

This sensor is suitable for measurements of low conductivities (0 $\mu\text{S}/\text{cm}$ to 300 $\mu\text{S}/\text{cm}$), e.g., in deion. water or for measurements in accordance with USP 645>.

The sensor can be connected to an 856 Conductivity Module with the adapter box (6.2103.160).

Phạm vi giao hàng 6.0914.040

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS

6.1236.050

Sleeve with SGJ 14/12 mm

Sleeve with SGJ 14/12 mm, polyethylene.



Phụ kiện tùy chọn

Order no.	Description	
6.1420.100	Flow-through measuring vessel For 12 mm dia. electrodes (pH, metal, conductivity)	
6.2324.110	Conductivity standard 100 µS/cm, 5 x 30 mL Conductivity standard for calibration of conductivity measuring cells with cell constant = 0.1/cm.	
6.2324.010	Conductivity standard 100 µS/cm 250 mL Conductivity standard for calibration of conductivity measuring cells with cell constant = 0.1/cm.	