



Hamilton PRP-X100 - 250/4.0

6.1005.010

The Hamilton PRP-X100 - 250/4.0 IC anion column is a robust separation column based on a polystyrene-divinylbenzene copolymer. It is mainly used with difficult matrices, e.g., dyes.

Phụ kiện tùy chọn

Order no.	Description
6.02821.000	Holder for the PRP-X100 Guard cartridge

Holder for the Hamilton PRP-X 100 Guard cartridge for connection to the analytical separation column.



6.1005.020

PRP-X100 guard column cartridge/4.0

For protecting Hamilton PRP-X100 columns. The cartridge efficiently removes contamination in the form of particles due, for example, to the growth of bacteria or algae.

