



Metrosep RP 3 Guard HC/4,0

6.1011.040

The Metrosep RP 3 Guard HC/4.0 is a guard column for universal use. It reliably protects the analytical separation column against contamination, securely removing lipophilic organic contamination, the smallest particles, traces of iron oxide and bacteria. The guard column is based on a polymer material and thanks to its larger pack volume, has a significantly higher capacity than the Metrosep RP 2Guard/3.5. It can be used throughout the entire pH range.

Phạm vi giao hàng 6.1011.040

Qt.	Order no.	Description
1 PCS	6.2744.060	Threaded stopper
		For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.

