



# Metrosep A Supp 16 Guard/2.0

6.1031.600

The Metrosep A Supp 16 Guard/2.0 reliably protects the Metrosep A Supp 16 analytical separation columns with 2 mm inner diameter against contamination. Thanks to the "On Column Guard System", the guard column is very easy to handle. The guard column screws easily onto the analytical column. No tools are required.

## Phạm vi giao hàng 6.1031.600

| Qt.   | Order no.  | Description      |
|-------|------------|------------------|
| 1 PCS | 6.2744.060 | Threaded stopper |

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.

