




Pt1000 temperature sensor for 909 UV Digester

6.1110.010

Pt1000 temperature sensor (class B) made of glass with fixed cable (0.5 m) for connecting to a 909 UV Digester.

Phạm vi giao hàng 6.1110.010

Qt.	Order no.	Description
1 PCS	6.1236.020	Sleeve with SGJ 14/12 mm
		Sleeve with SGJ 14/12 mm and O-ring.



Together with 6.2008.050 storage vessel holder, it provides a support for the electrode on 807 Dosing Units.

