



Pt1000 temperature sensor for CVS (installation length 9 cm)

6.1110.120

Pt1000 temperature sensor (class B) for Professional CVS Instruments.

Phạm vi giao hàng 6.1110.120

Qt.	Order no.	Description
100 Stk	6.1236.020	Sleeve with SGJ 14/12 mm
		Sleeve with SGJ 14/12 mm and O-ring.



Phụ kiện tùy chọn

Order no.	Description
6.2104.140	Electrode cable / 1 m / 2 x 2 mm For connecting temperature sensors with Metrohm plug-in head G with 780, 781, 826, 827 pH/Ion meters and all Titrandos

