



Pt100 temperature sensor for stability measuring devices (17.5 cm)

6.1111.010

Pt100 temperature sensor (length 17.5 cm) with fixed cable and mini-DIN-plug. This glass-free temperature sensor is well suited for use with a 711 Liquino and instruments for stability measurement (743 Rancimat, 763 PVC Thermomat, 892 Professional Rancimat, or 895 Professional PVC Thermomat).

Phụ kiện tùy chọn

Order no.	Description
6.2042.040	Clamp for temperature sensor
	Holding clamp for Pt100 temperature sensor for stability measuring instruments

