



Pt100 temperature sensor for stability measuring devices (30 cm)

6.1111.020

Pt100 temperature sensor (length 30.0 cm) with fixed cable and mini-DIN-plug. This glass-free temperature sensor is well suited for use with instruments for stability measurement (873 Biodiesel Rancimat and 893 Professional Biodiesel Rancimat).

Phụ kiện tùy chọn

Order no.	Description
6.2042.040	Clamp for temperature sensor
	Holding clamp for Pt100 temperature sensor for stability measuring instruments

