



Pt1000 temperature sensor made of stainless steel

6.1114.010

Pt1000 temperature sensor (class B) with fixed cable (1.2 m), with a needle made of stainless steel and a PEEK shaft. This glass-free temperature sensor is well suited for measurements:

- in semisolid materials (e.g., cheese, meat)
- in non-oxidizing materials at pH 1 to 13

Phạm vi giao hàng 6.1114.010

Qt.	Order no.	Description
1 PCS	6.1236.020	Sleeve with SGJ 14/12 mm
		Sleeve with SGJ 14/12 mm and O-ring.

