



1 mm platinum electrode tip for CVS

6.1204.190

Metal electrode tip made of platinum (Pt), diameter of electrode disk 1 mm ± 0.02 mm, polished surface, shaft made of glass, M3 thread connection. Concentricity error: 0.2 mm. Temperature range: 0 - 50 °C.

Phụ kiện tùy chọn

Order no.	Description
6.2103.110	Contact pin M3, 2 mm For 6.1241.060 Electrode shaft and 6.1204.1XX Electrode tips. With M3 thread. Used with VA equipment sets and for ELCD applications.

