



Glassy carbon electrode tip

6.1204.300

Electrode tip made of glassy carbon (GC) for Autolab RDE, electrode disk diameter 3 mm, M4 thread connection. Temperature range: 0 - 40°C.

Phạm vi giao hàng 6.1204.300

Qt.	Order no.	Description
1 PCS	6.2802.000	Polishing set for solid-state electrodes
		Polishing set with 1 polishing cloth and about 2 g aluminum oxide powder (grain size 0.3 µm)



Phụ kiện tùy chọn

Order no.	Description	
6.2103.120	Contact pin M4 / 2 mm For 6.1241.060 Electrode shaft and 6.1204.3XX Electrode tips. With M4 thread. Used with VA equipment sets.	