



Electrode holder

6.1241.020

Electrode holder for VA auxiliary electrodes with lock bolt. For use with glassy carbon rod (6.1247.000) or platinum pin (6.1247.010).

Phụ kiện tùy chọn

Order no.	Description
6.1247.000	Glassy carbon rod (GC) For electrode holder 6.1241.020 or 6.1241.120 for VA auxiliary electrode.



6.1247.010

Platinum pin (Pt)

For electrode holder 6.1241.020 or 6.1241.120 for VA auxiliary electrode.

