



Electrode holder

6.1241.120

Electrode holder for VA auxiliary electrodes with plug contact 5.5. mm. For use with glassy carbon rod (6.1247.000) or platinum pin (6.1247.010).

Phụ kiện tùy chọn

Order no.	Description
6.1247.000	Glassy carbon rod (GC)

For electrode holder 6.1241.020 or 6.1241.120 for VA auxiliary electrode.



6.1247.010

Platinum pin (Pt)

For electrode holder 6.1241.020 or 6.1241.120 for VA auxiliary electrode.



6.2112.130

Adapter plug contact 5.5 mm / screw contact, gold-plated

For VA electrodes

