





MME measuring head for Professional VA instruments

6.1256.020

Measuring head for operation with the Multi-Mode electrode pro or the scTRACE Gold.

Phạm vi giao hàng 6.1256.020

Qt.	Order no.	Description
1 PCS	6.2069.000	Park station for the measuring head of Professional VA/CVS instruments






Set with 1 stopper for electrode opening (e.g. for pH electrode or Pt1000 temperature sensor), 6 stoppers for openings (auxiliary solutions and/or automatic sample addition).



10 pcs, for Professional VA instruments, including wrench.



Phụ kiện tùy chọn

Order no.	Description	
6.1454.130	Form seal made of FPM for 61414xxx for KF and VA instruments	

