





SPE measuring head for Professional VA instruments

6.1256.030

Measuring head for operation with *screen-printed electrodes* (SPE) or the scTRACE Gold.

Phạm vi giao hàng 6.1256.030

Qt.	Order no.	Description
1 PCS	6.2069.000	Park station for the measuring head of Professional VA/ CVS instruments






Set with 1 stopper for electrode opening (e.g. for pH electrode or Pt1000 temperature sensor), 6 stoppers for openings (auxiliary solutions and/or automatic sample addition).



10 pcs, for Professional VA instruments, including wrench.



Phụ kiện tùy chọn

Order no.	Description	
6.1454.130	Form seal made of FPM for 61414xxx for KF and VA instruments	
6.2709.130	Adapter B-14/15 - SPE measuring head insert (VA) Adapter for use of the electrode shaft scTRACE Gold 6.1241.080 or the bismuth drop electrodes 6.0346.000 in the 884 Professional VA with SPE measuring head 6.1256.030 or in the 797 VA Computrace or 663 VA Stand with measuring head insert 6.1412.020.	