




RDE measuring head for Professional VA/CSV instruments with high chemical resistance


6.1256.040

Measuring head for operation with rotating disk electrodes. Solvent-resistant version for applications in nonaqueous media.

Phạm vi giao hàng 6.1256.040

Qt.	Order no.	Description
1 PCS	6.1244.050	Drive belts for rotating disk electrodes

For Professional VA/CSV instruments (3 pcs.).



1 PCS


6.1829.030

PTFE tubing / M6 / 21 cm, tapered



1 PCS

6.2069.000

Park station for the measuring head of Professional
VA/CVS instruments

1 PCS

6.2709.110

Stopper for the electrode opening on Professional VA
/CVS instruments

Set with 1 stopper for electrode opening (e.g. for pH electrode or Pt1000 temperature sensor), 6 stoppers for openings (auxiliary solutions and/or automatic sample addition).



Phụ kiện tùy chọn

Order no.	Description	
6.1454.130	Form seal made of FPM for 61414xxx for KF and VA instruments	

