



KF titration vessel lid with 5 openings

6.1414.030

For volumetric KF titrations.

Phạm vi giao hàng 6.1414.030

Qt.	Order no.	Description
1 PCS	6.1454.130	Form seal made of FPM for 61414xxx
		for KF and VA instruments



Phụ kiện tùy chọn

Order no.	Description	
6.1244.040	Set of seals Set of seals for 6.1414.030. The set is comprised of one seal and various O-rings.	
6.2730.010	Screw nipple Made of PP. For 6.1448.010 Septum or for mounting 6.1403.040 Drying tube or 6.0338.100 Electrode in 6.1414.030 and 6.1414.040 KF Titration vessel lids. Used for KF titrations.	
6.2730.020	Septum stopper With E.301.0041 O-ring. For 6.1448.010 Septum and 6.1414.030 Titration vessel lid. Used for KF applications.	