



Reaction vessel long for stability measurements

6.1429.050

Reaction vessel for stability measurements involving volatile oils, mineral oils and biodiesel blends with the Biodiesel Rancimat. In combination with 6.2418.130. Set of 100 pieces.

Phụ kiện tùy chọn

Order no.	Description	
6.1808.090	Thread adapter M8 outer / M6 inner Outer thread M8, inner thread M6, e.g. for mounting 6.2418.XXX Air tubes in 6.1429.XXX Reaction vessels of stability measuring instruments.	

6.2418.130

Air tube long for biodiesel measurements

Set of 100 pieces. In combination with 6.1429.050.

