



PEEK SGJ 29 stopper for KFC titration cell

6.1446.240

Stopper for volumetric titration in a titration cell (coulometry)

Phạm vi giao hàng 6.1446.240

Qt.	Order no.	Description
1 PCS	6.1543.200	Titration tip / M6 thread

Titration tip with antidiffusion valve.




1 PCS

6.01406.010

Adsorber tube with lid

Adsorber tube with lid for OMNIS Karl Fischer titrator, volumetric



1 PCS

6.02709.010

Stopper

Plastic stoppers, suitable for 2.1001.4220 OMNIS Karl Fischer Titrator.



1 PCS

6.02709.030

Screw nipple

Plastic screw nipple, suitable for 6014050X0 titration vessel lids.



1 PCS

6.2713.010

PTFE sleeve SGJ29

Sleeve for fat-free ground-joint connection with the standard ground-joint SGJ 29.

