




Titration head / 2 x SGJ 14

6.1458.070

Titration head for 75 mL beakers. With 2 openings SGJ 14, 2 openings 6.4 mm, and 3 pre-mounted dosing tips. Suitable for rinsing and aspiration equipment.

Phạm vi giao hàng 6.1458.070

Qt.	Order no.	Description
4 PCS	6.1808.000	Connection bushing M6



3 PCS

6.2726.090

Antidiffusion valve

For antidiffusion buret tips 6.1541.XXX, 6.1543.130, 6.1543.140,
6.1543.150, 6.1543.160, 6.1543.200



2 PCS

6.2740.020

Spray nozzle

For the fine-spraying of the rinsing solution.



Phụ kiện tùy chọn

Order no.	Description	
6.1829.050	Replacement tubings for retracted tubings in titration heads and robotic arms Three replacement tubings with tapering tip.	