



Robotic arm with titration head for 786 Swing Head, left or right swinging

6.1462.050

Titration head for Swing Head

Phạm vi giao hàng 6.1462.050

Qt.	Order no.	Description
2 PCS	6.2726.090	Antidiffusion valve
		For antidiffusion buret tips 6.1541.XXX, 6.1543.130, 6.1543.140, 6.1543.150, 6.1543.160, 6.1543.200





Phụ kiện tùy chọn

Order no.	Description	
2.802.0020	802 Stirrer for sample changer Rod stirrer for sample changer and Sample Processor. With 104 mm 6.1909.020 stirring propeller and permanently installed cable.	
6.1829.050	Replacement tubings for retracted tubings in titration heads and robotic arms Three replacement tubings with tapering tip.	