



Robotic arm with titration head and beaker sensor for 786 Swing Head, left or right-swinging

6.1462.150

Titration head with beaker sensor for 786 Swing Head, 2 x M10, 1 x SGJ 14.5, 2 x buret tip

Phạm vi giao hàng 6.1462.150

Qt.	Order no.	Description
3 PCS	6.1808.000	Connection bushing M6



2 PCS

6.2726.090

Antidiffusion valve

For antidiffusion buret tips 6.1541.XXX, 6.1543.130, 6.1543.140,
6.1543.150, 6.1543.160, 6.1543.200



Phụ kiện tùy chọn

Order no.	Description	
2.802.0020	802 Stirrer for sample changer Rod stirrer for sample changer and Sample Processor. With 104 mm 6.1909.020 stirring propeller and permanently installed cable.	
6.1829.050	Replacement tubings for retracted tubings in titration heads and robotic arms Three replacement tubings with tapering tip.	