



807 Dosing Unit ETFE 2 mL

6.1575.120

807 Dosing Unit ETFE 2 mL for dosing alkali reagents or hydrofluoric acid (HF) with Dosinos. Less precision and reproducibility is to be expected in comparison with glass cylinders.

Phạm vi giao hàng 6.1575.120

Qt.	Order no.	Description
1 PCS	6.1446.040	Threaded stopper / M6



1 PCS


6.2803.010

Grease (2 g)

Special quality (without silicone). For grease-free ground joint connections see 6.2713.XXX



Phụ kiện tùy chọn

Order no.	Description	
6.1574.120	Cylinder element glass 2 mL Glass cylinder element complete with cylinder, piston, cylinder base, valve disk, and centering tube. For 807 Dosing Unit 2 mL.	
6.1566.120	Cylinder element ETFE 2 mL For dosing units used for dosing alkali reagents or hydrofluoric acid (HF) with Dosinos. Less precision and reproducibility is to be expected in comparison with glass cylinders.	