



807 Dosing Unit ETFE 10 mL

6.1575.210

807 Dosing Unit ETFE 10 mL for dosing alkali reagents or hydrofluoric acid (HF) with Dosinos. Less precision and reproducibility is to be expected in comparison with glass cylinders.

Phạm vi giao hàng 6.1575.210

Qt.	Order no.	Description
1 PCS	6.1446.040	Threaded stopper / M6



1 PCS

6.2803.010

Grease (2 g)

Special quality (without silicone). For grease-free ground joint connections see 6.2713.XXX



Phụ kiện tùy chọn

Order no.	Description	
6.1574.210	Cylinder element glass 10 mL Glass cylinder element complete with cylinder, piston, cylinder base, valve disk, and centering tube. For 807 Dosing Unit 10 mL.	
6.1566.210	Cylinder element ETFE 10 mL For dosing units used for dosing alkali reagents or hydrofluoric acid (HF) with Dosinos. Less precision and reproducibility is to be expected in comparison with glass cylinders.	