



807 Dosing Unit ETFE 20 mL

6.1575.220

807 Dosing Unit ETFE 20 mL for dosing alkali reagents or hydrofluoric acid (HF) with Dosinos. Less precision and reproducibility is to be expected in comparison with glass cylinders.

Phạm vi giao hàng 6.1575.220

Qt.	Order no.	Description
1 PCS	6.1446.040	Threaded stopper / M6



1 PCS

6.2803.010

Grease (2 g)

Special quality (without silicone). For grease-free ground joint connections see 6.2713.XXX



Phụ kiện tùy chọn

Order no.	Description	
6.1574.220	Cylinder element glass 20 mL Glass cylinder element complete with cylinder, piston, cylinder base, valve disk, and centering tube. For 807 Dosing Unit 20 mL.	
6.1566.220	Cylinder element ETFE 20 mL For dosing units used for dosing alkali reagents or hydrofluoric acid (HF) with Dosinos. Less precision and reproducibility is to be expected in comparison with glass cylinders.	