



807 Dosing Unit ETFE 50 mL

6.1575.250

807 Dosing Unit ETFE 50 mL for dosing alkali reagents or hydrofluoric acid (HF) with Dosinos. Less precision and reproducibility is to be expected in comparison with glass cylinders.

Phạm vi giao hàng 6.1575.250

Qt.	Order no.	Description
1 PCS	6.1446.040	Threaded stopper / M6



1 PCS

6.2803.010

Grease (2 g)

Special quality (without silicone). For grease-free ground joint connections see 6.2713.XXX



Phụ kiện tùy chọn

Order no.	Description	
6.1574.250	Cylinder element glass 50 mL Glass cylinder element complete with cylinder, piston, cylinder base, valve disk, and centering tube. For 807 Dosing Unit 50 mL.	
6.1566.250	Cylinder element ETFE 50 mL For dosing units used for dosing alkali reagents or hydrofluoric acid (HF) with Dosinos. Less precision and reproducibility is to be expected in comparison with glass cylinders.	