



## 807 Dosing Unit 20 mL without accessories

6.1580.220

807 Dosing Unit with integrated data chip with 20 mL glass cylinder and light protection, mountable to a reagent bottle with ISO/DIN GL 45 glass thread.

Phạm vi giao hàng 6.1580.220

| Qt.   | Order no.  | Description                                                                                 |
|-------|------------|---------------------------------------------------------------------------------------------|
| 1 PCS | 6.2803.010 | Grease (2 g)                                                                                |
|       |            | Special quality (without silicone). For grease-free ground joint connections see 6.2713.XXX |

