




Bottle attachment / GL 45 / complete

6.1602.125

Enables direct dosing from reagent bottles with standardized GL 45 thread.

Phạm vi giao hàng 6.1602.125

Qt.	Order no.	Description
1 PCS	6.1602.105	Bottle attachment / GL 45 / green

For direct dosing from reagent bottles with GL 45 thread.



1 PCS

6.1446.080

Threaded stopper / M8



1 PCS

6.1805.080

FEP tubing / M6 / 25 cm

With light and kink protection.



1 PCS

6.1819.020

Tube

For M6 threaded holes. For 6.1602.1XX bottle attachments



1 PCS

6.2023.020

Clip for SGJ 14/15

Clip for SGJ 14/16

