




Bottle attachment / S 40 / complete

6.1602.135

Enables direct dosing from reagent bottles with S 40 thread (Merck...).

Phạm vi giao hàng 6.1602.135

Qt.	Order no.	Description
1 PCS	6.1446.080	Threaded stopper / M8



1 PCS

6.1602.115

Bottle attachment / S 40 / green

For direct dosing from reagent bottles with S 40 thread (Merck...).



1 PCS

6.1805.080

FEP tubing / M6 / 25 cm

With light and kink protection.



1 PCS

6.1819.020

Tube

For M6 threaded holes. For 6.1602.1XX bottle attachments



1 PCS

6.2023.020

Clip for SGJ 14/15

Clip for SGJ 14/16

