



Bottle attachment GL 45 for eluents and reagents

6.1602.160

For eluent, reagent, and rinsing bottles, with connectors for the adsorber tube and the aspiration tubing.

Phạm vi giao hàng 6.1602.160

Qt.	Order no.	Description
1 PCS	6.1446.040	Threaded stopper / M6



1 PCS

6.1602.105

Bottle attachment / GL 45 / green

For direct dosing from reagent bottles with GL 45 thread.



1 PCS

6.1446.080

Threaded stopper / M8

