



Sample rack, 12 x 6 mL, (Al)


6.2041.210

Sample rack for 889 IC Sample Center, suitable for 12 sample vessels. The following sample vessels can be used: 6.2419.000.

Material: Aluminum (Al)

Phụ kiện tùy chọn

Order no.	Description
6.1448.050	Aluminum septum caps / 1000 pieces Including Al sealing.



6.2419.000 Sample glass, 6 mL, 1000 pieces

Sample glass for up to 6 mL sample, suitable for 6.1448.050 septum, 1,000 pieces



6.2621.110 Crimping tongs

Crimping tongs for sealing 6 mL sample vials (6.2419.000, 6.2419.007).



6.2419.007 Sample glass, 6 mL, 100 pieces

Sample glass for up to 6 mL sample, suitable for 6.1448.057 septum, 100 pieces



6.2420.007 Sample glass with thread, 8 mL, 100 pieces

Sample glass with thread, 8 mL, suitable for 6.1448.067 screw cap, 100 pieces



6.1448.067

Screw cap with septum, 100 pcs.

Screw cap with septum, suitable for sample glass 8 mL 6.2420.007, 100 pieces



6.1448.057

Aluminum septum caps / 100 pieces

Including Al sealing.

