



Sample rack, 14 x 200 mL, (PVC)

6.2041.370

Sample rack for 810 Sample Processor, 814 USB Sample Processor, 855 Robotic Titrosampler "Light", suitable for 14 sample beakers. The following sample beakers can be used: 6.1459.310.

Plastic: Polyvinyl chloride (PVC)

Phụ kiện tùy chọn

Order no.	Description
6.1459.310	Sample beaker, disposable, plastic (PP), 200 mL, 1000 pieces
	Sample beaker, disposable, plastic (polypropylene, PP), 200 mL, suitable for 6.02041.020 OMNIS Software sample rack, 6.2041.370 und 6.2041.830 sample racks, 1000 pieces.

